## Arm Type Blood Pressure Monitor

Model: # BPM67

### INSTRUCTION MANUAL

PLEASE READ THIS INSTRUCTION MANUAL COMPLETELY REFORE OPERATING THIS UNIT

### FN













(♣) CE X X MD X IP 22 RoHS REACH



MDSS GmbH Schiffgraben 41 30175 Hannover Germany



**AViTA Corporation** 9F., No. 78, Sec. 1, Kwang Fu Rd., San Chung Dist., New Taipei City 24158, Taiwan, China



China, AViTA(WUJIANG)



Importer



Distributor

### Contents

Intended Use	1
Important Information before Use	2
Product Identification	6
Description of LCD Display	7
Battery Installation	8
Setting the Date and Time	9
Placement of the Pressure Sleeve	10
Measurement of Pulse Rate and Blood Pressure	12
Hypertension Indicator	13
Irregular Heartbeat Detector	14
Atrial Fibrillation (AFIB)	14
Memory Function	16
Error Codes	18
Troubleshooting	19
Cleaning and Disinfecting	20
Technical Specification	21
EMC Tables	22
Explanation of Symbols	24

#### Intended Use

The product automatically measures human being's Systolic, Diastolic blood pressure and pulse rate by oscillometric method. The measurement results are displayed on the LCD. Measurement position is at human being's upper arm. The intended use of this over-the-counter device is for adults with upper arm circumference ranging from 220 mm to 420 mm (Approx.8.7 ~ 16.5 inches) and for home use. When the device detects the appearance of irregular heartbeats such as atrial or ventricular premature beats during measurement, an indicated symbol will appear with measuring readings.

This device is designed only for adults.

### Type of Use/ Reuse

Multiple patient multiple use

#### Intended User

The Arm type blood pressure monitor is intended or both professional and consumer, and the patient is the intended operator. Patient selection criteria: Handicapped persons and children are the exception, as of handicapped persons and children need assistance by another person to use the device.

#### Contra-indications

Do not use in these cases: common arrhythmias such as atrial or ventricular premature beats, atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, age, renal diseases, patient motion, trembling, shivering.

### Important Information before Use

- Blood pressure measurements should only be interpreted by a
  physician or a trained health care professional who is familiar with
  your medical history. Through regular use of this device and
  recording of your measurements, you can keep your physician
  informed of the changes in your blood pressure.
- Perform your measurement in a quiet place. You should be seated in a relaxed position.
- Avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to making a measurement. If you are exhibiting signs of stress, avoid taking your measurement until the feeling subsides.
- 4. Rest 15 minutes prior to taking a reading.
- 5. Remove any constrictive clothing or jewelry that may interfere with the cuff placement.
- 6. Keep the monitor stable during measurement to achieve an accurate reading. Remain still; do not talk during the measurement.
- 7. Record your daily blood pressure and pulse readings on a chart.
- 8. Take your readings at the same time each day or as recommended by your physician to get an accurate indication of change in your true blood pressure.
- Wait a minimum of 15 minutes between readings to allow blood vessels to return to normal. The wait time may vary depending on your individual physiological characteristics.
- 10. Although such cases are rare, for those with an extremely weak pulse or irregular pulse, errors may result which prevent proper measurement. If abnormal variations are noticed, consult with your physician or trained healthcare professional.
- 11. This device is intended for adult use. While taking a measurement, you can stop the inflation or deflation process of the cuff at any time by pressing the POWER button.

- 12. If any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.
- 13. Do not expose the device to extreme temperatures, humidity dust or direct sunlight as this may cause it to malfunction.
- 14. Please comply with the storage and operating conditions defined in 'Technical Specification'. Storing or using the device outside of the specified temperature and humidity range can affect measurement accuracy or the function of the device.
- 15. If the device was not stored within the minimum/maximum permissible storage conditions, a waiting period of at least 2 hours must be observed before using it under the specified operating conditions ( 'Technical Specification' ) or an ambient temperature of approx. 20 °C.

For Customer Service, the blood pressure monitor is calibrated when manufactured; it is recommended that the accuracy should be maintained and calibrated by manufacture triennially (every 3 years). To obtain the service please contact AViTA Corp. for the address of the repair location. Enclose with the Proof of Purchase. Include \$10.00 USD for the return shipping and handling. Accompany with a letter, with your name, address, phone number, and description of the specific problem or routine check-up. Pack the device carefully with bubble wraps (if there is) to prevent damages cause during transit. Due to possible losses in transit, it is recommended insuring the device with return receipt requested. If in any way of assistance of setting up, using, maintaining or to report unexpected operation/adverse events please contact manufacturer or local representative for further assistance.

## A CAUTION:

- Beware of blood flow interference and resulting harmful injury to the patient caused by continuous cuff pressure.
- While using monitoring ME equipment and blood pressure monitor simultaneously being used on the same limb, pressurization of the cuff may cause monitoring ME equipment temporarily dysfunction.
- Reading can be affected by the measurement site, the position of the patient, exercise, or patient's physiologic condition.
- Automated sphygmomanometer can be affected by extremes of temperature, humidity and altitude.
- Frequent measurements can cause injury to the patient due to blood flow interference.
- Improper operation of automated sphygmomanometer may result in prolonged impairment of patient blood circulation, please read manual for proper use, if still not clear, contact professional health care personnel or local distributor for assistant.
- · Retake the measurement if unexpected readings are obtained.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g., intravascular access or therapy, or an arteriovenous (AV) shunt.
- This product is suitable for use in the home healthcare environment.
- Keep this device out of the reach of children. Strangulation may result from baby or child entanglement in cables.
- Please keep this device away from pets, pests, and children.
- Preventing potential allergic reaction, please avoid the device in direct contact to patient's wound.

- Do not use the cuff on people who have undergone a mastectomy.
- Do not apply the cuff over a wound, as this can cause further injury.
- Do not apply the cuff other then the original manufacturer provided.
- Do not use in these cases (e.g. Device for use in an ambulance, helicopter or professional environment)
- Cuff pressure 0 300 mmHg
- Reduction rate: ≤30S Refer to IEC 80601-2-30
- No modification of this equipment is allowed.
- · High BMI health condition users may result varied.
- Users shall notify serious adverse event to the central competent authority or its commissioned agency, legal entity, or manufacture.

### **Product Identification**



- 1 Display
- 2 Cuff
- 3 Memory button (M key)
- 4 Start/Stop button (Power key)
- 5 Set button (S key)

## **Description of LCD Display**

Low battery indicator

WHO indicator

88÷88

Date & Time

Irregular Heartbeat Symbol

**A** 

AF Symbol

M

Memory Symbol

188

Memory Set

V

Heartbeat Symbol

188

Pulse rate

**\$** 

Inflation / Releasing air

Systolic Pressure

888

Diastolic Pressure

12

User Symbol


### **Battery Installation**

### Low battery warning:

It is necessary to replace the batteries when the Low Battery symbol "[]" appears on the display, or when the display does not turn on after the POWER key is pressed.

### Replacing the Battery:

- 1. Press down on latch and lift the cover on the bottom of the monitor.
- Insert or replace 4x 1.5 V AAA batteries into the battery compartment, ensuring to match the indicated polarity symbols. Always use new batteries.
- 3. Replace the battery cover.



# **↑** CAUTION:

Risk of explosion if battery is replaced by an incorrect type. Dispose of used batteries according to the instructions.

### NOTE: Battery-operated

- Please properly dispose of the batteries away from small children and heat. Avoid the children accidentally swallow the battery.
- 2. It is recommended to remove the batteries if the device will not be used for more than 1 month.
- 3. Batteries must be disposed of in accordance with local environmental and institutional policies.
- 4. It is recommended not to use rechargeable, unqualified or different spec battery may damage the device or cause circuit shortcut.

### Setting the Date and Time

It is necessary to set the date and time for the device every time batteries are initially installed or replaced.

To set the date and time, proceed as follows:

- While in power off mode, press and hold the "S" key for at least 3 seconds to enter Date and Time setting procedure and the Year value will begin to flash.
- Press the "S" key to advance the display to the desired year, press the "M" key to confirm the year.
- Next, the month value will blink. Repeat step 2 to set the month and date, then hours, then minutes.
- After setting the minutes, the device will automatically exit out of the date/time setting mode and shut off.

#### Placement of the Pressure Sleeve

It is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from

- Remove any constrictive clothing or jewelry that may interfere with cuff placement.
- 2. Be seated at a table or desk with your feet flat on the floor.

#### Note:

Blood pressure naturally varies from one arm to the other; therefore, measure your blood pressure on the same arm to ensure comparability of the two readings.

- 3. Fully open the cuff by pulling the bottom of it. Then grip the clip of the device to widen the cuff, creating a clamp shape.
- 4. Insert your left arm into the created clamp shape, with the direction indicator align with your brachial artery.
- 5. The bottom edge of the cuff should be positioned approximately one inch above the elbow joint.
- Reach underneath your left arm with your right hand, and pull the end of the cuff towards your body to tighten the cuff. Wrap and secure the cuff.
- The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger easily between your arm and the cuff.

### Note:

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.

#### Measurement of Pulse Rate and Blood Pressure

- Press the Power key to turn the power on. After full display is shown, the values for the last reading will appear on the display. If there is no measurement, the unit displays the value "0".
- 2. AF mode selection will appear after the value of last reading. Press the Power key again to start the measurement.
- After the self-test, the blood pressure monitor starts to measure. The cuff will automatically begin to inflate, with the display showing the increasing pressure in the cuff.
- As the pressure increases, the indicator will increase upwards according to the pressure value on the display.
- 5. As the cuff inflates, the monitor automatically determines your ideal inflation level. This monitor detects your blood pressure and pulse rate during inflation.
- To detect the heartbeat, the heartbeat symbol will appear and continuous flashes on the LCD display.
- Your blood pressure measurement and pulse will display simultaneously on the screen.
- The Hypertension Indicator will indicate your reading range on the display separately.
- Press the Power key to turn the device off, conserving energy and battery life. The device will automatically shut-off after approximately 2 minutes of idling.


If there is an Error message after the measurement, please refer to Page 18 Error Codes Table.

#### Note:

- Users are recommended to use BPM67 daily at a fixed time, frequency of usage from 1 to 3 times can be suggested from physician or doctor's advice.
- Base on individual body condition the typical operation time takes approximately 30 seconds to 1 minute.

## Hypertension Indicator

This device features our unique Hypertension Indicator. The World Health Organization has established globally accepted standards for the assessment of high or low blood pressure readings. The below chart should be considered only as a guideline, always consult with your physician or health care professional to interpret your individual results.



### Irregular Heartbeat Detector

Your digital blood pressure monitor features an Irregular Heartbeat Detector. This feature allows users to accurately monitor blood pressure even if an irregular heartbeat should occur. When an irregular heartbeat is detected, the icon "e" will appear on the display.

#### Note:

Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat if this symbol appears frequently.

#### IMPORTANT INFORMATION:

This blood pressure monitor is not designed for use by people with arrhythmias nor for diagnosing or treating an arrhythmia problem. As a safeguard, we recommend that if you have arrhythmias such as atrial or ventricular premature beats and atrial fibrillation or any other special conditions you should check with your physician before using your blood pressure monitor.

### Atrial Fibrillation (AFIB)

This device is a blood pressure monitor that also analyses heart rate variability during measurement.

While in power off mode, press Power key into select AF mode, and use "S" key to select "ON" or "OFF" the AF mode. Confirm your selection with the Power key to turn the power on and start to measure your blood pressure.

This feature allows users to accurately monitor blood pressure even if atrial fibrillation should occur. When atrial fibrillation is detected, the "" icon will appear on the display.

#### Note:

- Please consult with your physician or trained healthcare professional for further information regarding an atrial fibrillation if this symbol appears frequently.
- 2. Even if the AF symbol does not appear, there is still a possibility of Afib.

### **Memory Function**

Recalling Measurements in Memory:

You can recall up to 100 measurements for each user, plus those average values stored measurements in memory to share with your physician or trained healthcare professional.

If you press the M key, the unit will first display the average of last 3 currently stored measurements. Continue to press the M key to successively view the average value of the morning measurements and evening measurements for the last 7 days will be displayed.

The duration of morning is AM5:00 – AM9:00, The duration of evening is PM6:00 – PM8:00

If you press the M key again, the measurements will appear on the display from most current to oldest. And the memory number will appear on the display.

All results for a given measurement will display, including measurement results, pulse rate, Hypertension Indicator, Irregular Heartbeat alert, and date/time stamp.

When the number of readings exceeds 100, the oldest data will be replaced with the new record.

Press the Power key to turn the device off at any time when you check the stored measurements

### Clearing Measurements from Memory:

From the memory mode, press and hold down the M key until the display shows CLr. This indicates that all measurements have been erased

### Two users with one device

This device is designed to be capable of recording the measurements of two people separately. To switch to the other user, please follow the steps below:

- 1. Press the "S" key while in power off mode and the user number will appear below the Date & Time on the display.
- 2. Change the user number by pressing the "S" key.
- After selecting the user, press the power key to confirm and proceed to the measuring steps.
- 4. For checking the measurement records, press the "M" key after selecting the user.

## **Error Codes**

Err Code	Meaning	Corrective Action
Err 00	No pulse or detect pulses not enough.	Take off heavy clothes and retry again.
Err 01	The cuff is not fastened correctly, cuff pressure leakage or inflation too low or Overpressure protection.	The Arm cuff is not fastened properly. Re-apply the cuff, and take a measurement again.
Err 02	Inaccurate reading	Rest a while, relax and retry again.
Err 03	Inflation or deflation fail during the measurement	The Arm cuff is not fastened properly. Re-apply the cuff, and take a measurement again.
Err	Memory error.	Take off batteries to reboot the device, and then take another measurement. If the Err code persists after reinstalling the batteries, please contact the supplier or manufacturer.
	Low batteries	Replace all batteries with new ones.

# Troubleshooting

Nothing appears in Ra		
	atteries are drained.	Replace all batteries with
the display even		new ones.
	attery are not correctly	Reinsert batteries in the
turned on. ali		correct position.
		Replace all batteries with
Symbol appears.		new ones.
	colder temperatures	Warm up the batteries, or
	atteries have weaker	use the device in a
		warmer setting.
		Use Alkaline batteries and
time is		replace all batteries at the
inconsistent.		same time with same
inconsistent.		brand batteries.
No reading after Ba		Replace all batteries with
measurement.		new ones.
		Adjust patient and Arm
		cuff to measure.
		Rest a while, relax and
		measure again.
da	aries tilloughout the	ineasure again.
Suspicious pulse Bo	odily movement during	Refrain from moving
	evice use.	during measurement.
		Do not take
	ter exercise or exposure	
		exercise or coming back
		from the outdoors.
Power switches off Sy		
automatically.		Push the power button again, and then begin
automatically.		measure again.
During massuring It	could be a normal	Relax, and try to take a
	essure is higher than	measure agáin.
l th	le initial pressure value, le device automatically	
	umps to a higher	
l bi	ressure by 40mmHg	
	ach time.	Classic that the Auron souff is
	ne Arm cuff is not	Check that the Arm cuff is
ras		fastened properly and
		lretake the measurement. I

## Cleaning and Disinfecting

- Only use a soft, damp cloth to clean the monitor. Please do not use thinner, alcohol, detergents or solvents.
- The cuff can be cleaned carefully using a slightly damp cloth and mild soap solution. Do not completely immerse the cuff in water.
- It is recommended to clean and disinfect the cuff regularly or after each use, especially when used by several users, to prevent infection.
   The cuff should be disinfected, particularly on the inside, by wiping with a disinfectant. Use a disinfectant that is compatible with the cuff materials, e.g. 75 % ethanol or isopropyl alcohol.
- Keep the monitor in the appropriate carriage to protect it from external influences. And store it at appropriate condition.

## **Technical Specification**

· Measuring range:

Blood Pressure: 40~255 mmHg Pulse Rate: 40~199 beats/min

· Calibration Accuracy:

Blood Pressure : ± 3 mmHg Pulse rate : ± 4% of reading

· Operating environment:

10°C~40°C 10°C √40°C

15% to 85% relative humidity (non-condensing)

700-1060 hPa ambient pressure

• Storage/ Transportation environment :

-20 to 50 °C ₂₂υς 1 50°C

15% to 85% relative humidity (non-condensing)

700-1060 hPa ambient pressure

• Power Source: 4 x 1.5 V AAA batteries

Body Weight:

approx. 250g (exclude batteries) +/- 5%

• Body Dimensions :

approx. L 130 mm x W 94 mm x H 33.6 mm

Cuff circumference (M/L Size):
 approx. 22 ~ 42 cm (9" ~ 17")

Lifetime: 3 years

• Expected service life: 10,000 measurements

Note: After 3 years life time or 10,000 measurements, device material may experience degradation, measurement accuracy may vary.

## **EMC Tables**

customer or	tne user of <b>BP</b> r	vi67 must make sure
		ctromagnetic emissions
facility 6	environment a)	ENVIRONMENT a)
CISPR 11		CISPR 11
		BLE Function)
		Group 2 Class B (With
	/	BLE Function)
		DC (t)
Only the AC input needs to be tested  Voltage fluctuations Not applicable		s to be tested
respiratory treatment facilities in hospital or clinics. The more restrictive		
acceptance limits of Group 1 Class B (CISPR 11) have been considered		
and applied. The equipment is suitable for use in the mentioned		
environments when directly connected to the Public Mains Network.		
licable in this	s environment	unless the ME
EQUIPMENT and ME SYSTEM used will be connected to the PUBLIC		
MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard.		
ırd.		
acturer's d	eclaration - Ele	ctromagnetic immunity
: FN4C	II	II-
on Basic EMC Immunity test levels		
		HOME HEALTHCARE ENVIRONMENT
tilletilou		EINVIROINIVIEINI
61000-4-2		<u> </u>
01000-4-2		•
	± ∠ KV. ±4KV	±. ±8 KV. ±15 KV all
61000-4-3		±, ±8 kV, ±15 kV air 10 V/m b) 80MHz - 2.7
	customer or an environma e	S Group 1 Class A (Not BL Function) Group 2 Class A (With BLE Function) Not applicable (Note: Power by Battery Only the AC input needs Not applicable (Note: Power by Battery Only the AC input needs in the AC input needs

Proximity fields from RF wireless communications equipment	IEC 61000-4-3	COMPLIANT NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the BPM67 can be requested from supplier using the contact information provided in this manual. However, it is advisable to keep the equipment at an adequate distance of, at least, 0.5 m from mobile phones or other RF communications transmitters to minimise possible interference.
RATED power frequency magnetic fields.	IEC 61000-4-8	30 A/m c) 50 Hz or 60 Hz

a) The equipment is suitable for use in Home Health Environments and Professional Health Care Environments limited to patient rooms and respiratory treatment facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits have been considered and applied.

b) Before modulation is applied.

c) This test level assumes a minimum distance of at least 15 cm between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic fields.

# **Explanation of Symbols**

Symbol	Definition
<b>C €</b> 2797	The CE marking with the Registration Number of the Notified Body. This denotes the compliance of Regulation (EU) 2017/745
MD	Medical Device
•••	Manufacturer
EC REP	Authorized representative in the European Community
~~ <u> </u>	Date of manufacture (YYYY-MM-DD or YYYY-MM)
LOT	Batch code (YYMMWWWW)
SN	Serial number (YYMWWWXXXXX)
<del>*</del>	Keep dry
*	Temperature limit

<u></u>	Humidity limitation
<b>♦•</b> ♦	Atmospheric pressure limitation
$\triangle$	Caution
<b>&amp;</b>	Consult the instruction for use
X	Disposal information: Should you wish to dispose of the article, do so in accordance with current regulations. Details are available from your local authority. WEEE 2012/19/EU Directives
RoHS	This product fulfilling the requirements of the RoHS Directive 2011/65/EU.
REACH	This product fulfilling the requirements of the REACH Directive EC 1907/2006 and its amendments, do not contain Substances of Very High Concern in concentration above the limit of 0.1 %. No substance(s) is/are present in the parts of the product above the concentration of 0.1 % weight by weight.
Ů	Stand-by
☀	Device classification type BF

IP 22	This product meets the basic safety and essential performance requirements indicated in the IP22 conditioning test (protection against solid foreign objects of 12.5mm Ø and greater and against vertically falling water drops when enclosure tiled up to 15°)
***	The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
	Importer
	Distributor
#	Model Number
₹ CN	Country of Manufacturer
UDI	Unique Device Identifier
紊	Keep away from sunlight

Electronic IFU available at http://www.avita.com.tw